

MOISTURE METER


DHS SERIES

DHS series electronic halogen moisture rapid tester adopts the principle of halogen lamp pyrolysis gravimetry and is a rapid moisture detection instrument. Halogen heating can reach the maximum heating power in the shortest time, and the sample is quickly dried at high temperature.



PRODUCT FEATURES

- 5-inch high-definition smart touch display, real-time dynamic analysis of curves, direct display of moisture percentage/dry weight percentage/dry weight;
- The ultimate heating temperature of 180°C meets the testing requirements of most samples;
- 9 sets of recent test records can be retrieved at any time and test results can be printed;
- Stainless steel heated inner cabin improves sample testing efficiency.

Product Model	LC-DHS-16A		LC-DHS-20A	
Image				
Voltage (V)	220			
Frequency (Hz)	50±1			
Power(W)	350			
Display Mode	5 "HD touch screen			
Max Weighing Value (g)	50		110	
Weighing Accuracy(g)	±0.002		±0.001	
Moisture Content Measurement Range (%)	0-100			
Temperature Range (°C)	50-180			
Temperature Control Accuracy (°C)	±1			
Timing Range (min)	0-99			
Working Size (mm)	Φ100			
Heating Method	High efficiency halogen lamp			