

WATER BATH


WB+ SERIES

WB+ series water bath adopts LCD and PID temperature control mode, with fast temperature rise and accurate temperature control. This series of thermostatic water bath pot has the function of preventing dry burning, which can avoid damage to the inner tank and heating pipe of the water bath pot caused by dry burning.



PRODUCT FEATURES

- With magnetic stirring function,it can accelerate the uniform water temperature;
- Connect the inner tank to the drainage pipe, and open the valve plug to automatically drain water (Limited to LCH-WB-4/6 models);
- LCD display, PID temperature control technology, fast temperature rise and accurate temperature control;
- It is high-temperature resistant, well sealed and durable;
- High quality temperature sensor is adopted, with accurate temperature measurement and long service life;
- Standard timing and anti dry burning function, safe and reliable.

Product Model	LCH-WB-8+
Image	
Voltage (V)	AC220
Frequency (Hz)	50
Power (W)	1800
Display Mode	LCD
Temperature Range (°C)	RT-100
Temperature Control Accuracy(°C)	±0.5
Timing Range	1-99h59min
Effective Volume (L)	30
Inner Dimension (mm)	685*390*225
Pipe Weep	√
Stirring Function	√