

# WATER BATH





## WB SERIES

WB series water bath adopts LCD and PID temperature control mode, with fast temperature rise and accurate temperature control. This series of thermostatic water bath pot has the function of preventing dry burning, which can avoid damage to the inner tank and heating pipe of the water bath pot caused by dry burning.



## PRODUCT FEATURES

- Connect the inner tank to the drainage pipe, and open the valve plug to automatically drain water (Limited to LCH-WB-4/6 models);
- The instrument adopts LCD display with clear and beautiful display contents, accurate and reliable, PID temperature control technology, fast temperature rise and accurate temperature control;
- The inner tank is made of high-quality stainless steel and stamped at one time. It is high-temperature resistant, well sealed and durable;
- High quality temperature sensor is adopted, with accurate temperature measurement and long service life;
- Standard timing and anti dry burning function, safe and reliable.

Product Model	LCH-WB-1	LCH-WB-2	LCH-WB-4	LCH-WB-6
Image				
Voltage (V)	AC110-220			AC220
Frequency (Hz)	50			
Power (W)	300	600	1000	1500
Display Mode	LCD			
Temperature Range (°C)	RT-100			
Temperature Control Accuracy(°C)	±1			
Timing Range	1-99h59min59s			
Effective Volume (L)	4.96	6.5	14	22
Inner Dimension (mm)	177*187*150	300*150*150	330*300*150	540*340*230
Pipe Weep	/	/	√	√