

# MINI VORTEX MIXER


## MX SERIES

MX series mini vortex mixer is characterized by its compact size and rapid and thorough mixing. The liquid is in a vortex shape and can mix all the test solutions attached to the tube wall. It is suitable for uniform mixing of liquid in general test tubes, centrifuge tubes, flasks, liquid separators, and funnels.



## PRODUCT FEATURES

- Small and compact appearance, easy to carry;
- Brushless DC motor, strong power and maintenance-free;
- Imported silicone vibrating head, strong toughness, sufficient flexibility and long service life;
- The sample can be swirled quickly within 5 seconds;
- Using flexible anti-vibration foot pads, metal base, and triangular structure, the instrument has good stability, low operating noise and no displacement during operation.

Product Model	LCH-MX-M2
Image	
Voltage	AC100-240V DC12V 1A
Frequency (Hz)	50-60
Power (W)	12
Speed Range (rpm)	500-3500
Circumferential Diameter (mm)	4.5
Max Throughput (ml)	50
Max Diameter of Test Tube (mm)	30
Press Vibration	Inching+Speed Regulation
Footing Material	Gel
Motor Type	DC Brushless Motor