

MULTI MODULE VORTEX MIXER


MX SERIES

MX series multi module vortex mixer is characterized by small size, mixing fast, thorough, liquid vortex can be attached to the wall of the test liquid all mixed evenly, suitable for general test tube, centrifugal tube, liquid, widely used in biochemistry, genetic engineering, medical and other experimental needs.



PRODUCT FEATURES

- LCD large screen display speed and time, easy to set read;
- Suitable for short time (point) or long time continuous work;
- Dc brushless motor, long life, strong power, maintenance-free special vibration head, strong toughness, elastic enough, long service life;
- Wide speed range, maximum speed 3000rpm, stepless speed regulation;
- The use of high quality pads, good stability of the instrument, low noise operation.

Product Model	LCH-MX-S				
Image					
Voltage	AC100-240V DC12V 2A				
Frequency (Hz)	50-60				
Power (W)	20				
Display Mode	Scale				
Speed Range (rpm)	200-3000				
Timing Range	/				
Circumferential Diameter (mm)	4.5				
Standard Head Diameter (mm)	25				
Press Vibration	Continuous Operation/Inching				
Motor Type	DC Brushless Motor				
Optional Modules	Φ4mm*48	Φ7.5mm*15	Φ9mm*16	Φ13mm*8	Φ17mm*8
Module Image	