

OVERHEAD STIRRER

OES SERIES

OES series overhead stirrer is suitable for mixing and stirring low, medium and high viscosity liquids in the laboratory. The stirring machine adopts high-quality engineering plastic as the body, the operating is stable, the output uses variable speed increase mechanism, so that the output torque is multiplied.

PRODUCT FEATURES

- The LCD screen displays the set speed and actual speed, which is intuitive and clear (LD items);
- DC permanent magnet motor, low noise, maintenance-free, stable operation;
- The overload will stop automatically to ensure safe operation;
- Stable start-up, effectively preventing sample overflow;
- Automatically memorize the final set working parameters, multi-segment programming control.



Product Model	LCH-OES-20LD
Image	
Voltage (V)	AC100-240
Frequency (Hz)	50 - 60
Power (W)	72
Display Mode	LCD
Max Stirring H ₂ O Volume (L)	20
Max Stirring Viscosity (mPa · s)	10000
Speed Range (rpm)	100-2000
Timing Range	1-99h59min
Max Torque (N · cm)	10
Motor Type	DC Permanent Magnet Motor